

# DIVIDED LEG SECTION DOUBLE JOINT SHORT

Divided leg section double joint, short with dual joint abduction for Promerix and new, smarter Practico operating tables.

Order no. 145480



<b>Bredd</b>	260 mm of one section
<b>Längd min-max</b>	555 - 655 mm
<b>Length (mattress base)</b>	625 mm
<b>Längd (table top)</b>	655 mm
<b>Bredd (table top)</b>	540 mm
<b>Viktintervall</b>	7,5 kg of one section
<b>Product Code</b>	145480