


# GRAND PROMERIX

Powerful and smooth Grand Promerix is heavy-duty, yet versatile operating table for all kind of surgical procedures. Improved ergonomic design improves the comfort of the whole surgical team. Comfort leads to concentration, confidence and better co-operation. That's efficiency.



## C-ARM ACCESS

The renewed slim base gives you perfect C-arm access. Grand Promerix meets even the strictest demands of modern operating theatres easily.



## IMPROVED ERGONOMICS

Grand Promerix has lots of improved details for better ergonomics, for example, hidden, smaller castors and inwards cutout for increased foot space. Human-centric design at its best.



### INCREASED HYGIENE

Details redesigned with patient safety in mind. Rethought frame make Grand Promerix hygienic and easy to clean.



### FLUENT LOCK™

Our Fluent locking system is both user-friendly and reliable. Connect with one click.



### HYGIENIC AND EASY TO CLEAN

New bellows and seat section make Grand Promerix hygienic and easy to clean.



### NEW DESIGN - IMPROVED DETAILS

Improved details improve hygiene and ergonomics for the operating room staff.

Titel	Sliding		Grand Promerix - Non-Sliding	
Längd (table top)	2100 mm		2100 mm	
Bredd (table top)	540 mm		540 mm	
Bredd (bas)	652 mm		652 mm	
Höjjustering	650 - 1110 mm		610 - 1070 mm	
Max. patientvikt	275 kg	600 lb	325 kg	700 lb
Säker arbetsbelastning (SWL)	275 kg	600 lb	325 kg	700 lb
Chockläge	-35 °		-35 °	
Anti-Chockläge	+35 °		+35 °	
Huvudstödsjustering	-45 - +45 °		-45 - +45 °	
Ryggsektionsjustering	-45 - +80 °		-45 - +80 °	
Bensektionsjustering	-105 - +70 °		-95 - +70 °	
Hjulstorlek	80		80 mm	
Primärspänning	100 - 240 V		100 - 240 V	
Sekundärspänning	24 V		24 V	
Max. effektförbrukning	700 W		700 W	
Frekvens	47-63 Hz		47-63 Hz	
Vikt	325 kg		300 kg	
Sidolutning	-25 - +25 °		-20 - +20 °	
Nedre ramlängd	1317 mm		1317 mm	
Max lyft kapacitet	490 kg		490 kg	
Product Code	100060000		100060000	