

K-TEK ARM BOARD

Integrated clamp snap locks onto rail at any position. Designed to prevent fluid ingress. Radiolucent, ultradurable strong material. Eliminates gear grinding and pinch hazards. One-hand adjustable. Pad not included.

Order no.

100011801 K-TEK2005 Arm Board, EU-rail

100011801U K-TEK2000 Arm Board, US-rail

100011801UK K-TEK2006 Arm Board, UK-rail

100012821 K-TEK2015 Arm Board Pad 3"/7,5 cm black



Max. patientvikt

227 kg

Vikt

1,8 kg

Product Code

100011801, 100011801U,
100011801UK, 100012821