

# K-TEK 250 STIRRUP PAIR 12814

The Kyra 250 is a junior stirrup designed for patients aged 6-12 years old but also effective in smaller/lighter adults weighing up to 250 lbs/113 kg.

- More working space with wider boot rotational range
- Lightweight and Strong with Kyra's Core-Grip™ technology

## Order information

100012814 K-TEK 250 STIRRUP PAIR

100012808 K-TEK™ STIRRUP CART

100012811 K-TEK™ STIRRUP BOOT PAD KIT

100 011 043 RAIL CLAMP, EUR and UK

100 011 043U RAIL CLAMP, US



<b>Max. patientvikt</b>	113 kg	250 lb
<b>Vikt</b>	4,6 kg	10,2 lb
<b>Movement angle</b>	+25	
<b>Lithotomy Range</b>	+85 - -55 °	
<b>Abduction range</b>	+25 - -9 °	
<b>Product Code</b>	100012814	