

K-TEK 350 STIRRUP PAIR I 2806

The K-Tek Comfort Stirrups for patient positioning for a wide range of gynecologic, urologic and laparoscopic procedures.

The Stirrups increase surgeons' operating room performance by permitting simple adjustment of lithotomy and abduction intra-operatively.

- Designed to relieve pressure with soft padding that encapsulates the foot.
- K-TEK™ silicone closure straps making cleaning easy.
- Flexible finned boots accomodates various leg sizes.
- Superior access with lithotomy range of motion from -55° to +85° aids in relieving clinican shoulder/arm strain.

Order information

I00012806 K-TEK COMFORT 350 STIRRUPS

I00012808 K-TEK™ STIRRUP CART

I00012811 K-TEK™ COMFORT™ 350 STIRRUP BOOT PAD KIT

I00 011 043 RAIL CLAMP, EUR and UK

I00 011 043U RAIL CLAMP, US





Max. patientvikt	159 kg	350 lb
Lithotomy Range	+85 - -55 °	
Abduction range	+35 - -9 °	
Product Code	100012806	