

K-TEK 500 STIRRUP PAIR 12812

The K-tek Comfort™ 500 boot accommodates heavier patients with smaller legs (patients weighing up to 500 lbs/227 kg).

Order information

I00012812 K-TEK 500 STIRRUP PAIR

I00012808 K-TEK™ STIRRUP CART

I00012811 K-TEK™ STIRRUP BOOT PAD KIT

I00 011 043 RAIL CLAMP, EUR and UK

I00 011 043U RAIL CLAMP, US



Max. patientvikt	227 kg	500 lb
Vikt	6,6 kg	14,6 lb
Movement angle	+25	
Lithotomy Range	+85 - -55 °	
Abduction range	+25 - -9 °	
Product Code	I00012812	