

K-TEK 800 STIRRUP PAIR 12813

The KYRA Comfort™ 800 boot accommodates heavier patients with smaller legs (patients weighing up to 800 lbs/363 kg)

Order information

I00012813 K-TEK 800 STIRRUP PAIR

I00012808 K-TEK™ STIRRUP CART

I00012811 K-TEK™ STIRRUP BOOT PAD KIT

I00 011 043 RAIL CLAMP, EUR and UK

I00 011 043U RAIL CLAMP, US



Max. patientvikt	363 kg	800 lb
Vikt	6,6 kg	14,6 lb
Movement angle	+/-25	
Lithotomy Range	+85 - -55 °	
Abduction range	+25 - -9 °	
Product Code	I00012813	