

# LIGHT NON-MOTORIZED LEG SECTION 60203

Light leg section improves work ergonomics of the OR staff. It is very light and easy to handle. Light weight leg section attaches to the operating table with easy locking mechanism.

Light leg section for non-motorized legs, 100060203, is suitable for New, smarter Practico operating table with non-motorized leg joints, Practico electro-mechanical and Practico Manual operating tables. To be used with Promerix operating table, 25- or 40-section is needed first to be attached to the seat section.

Includes 80 mm thick, antistatic VEF (viscoelastic foam) mattress.

Order code:

100060203 Light leg section for non-motorized legs, VEF mattress

For motorized leg joints, see Light leg section for motorized legs  
>>



Längd	630 mm
Bredd (table top)	540 mm
Höjd	76 mm
Säker arbetsbelastning (SWL)	40 kg
Vikt	4,9 kg
Product Code	100060203