

DORO HORSESHOE SYSTEM

Universal Horseshoe Headrest for adult or pediatric applications in prone or supine positioning. Mounted to a Doro Swivel Adaptor Aluminium or to the Doro Horseshoe Table Adaptor. Gel Pads interchangeable for adult or pediatric application. Permits vertical and lateral adjustments for proper patient positioning.



Horseshoe system

100 300 900 Adult headrest with Extension Bar
 100 300 800 Pediatric headrest with Extension Bar

Pads

100300910 Adult Right
 100300911 Adult Left
 100300810 Pediatric Right
 100300811 Pediatric Left

Cervical Spine Support

100301200

Product Code

100 300 900 100 300 800 100300910
 100300911 100300810 100300811
 100301200
