

# DORO TEFLON HEADREST SYSTEM

## Features:

- Autoclavable and machine washable system which conforms to WHO Infection Control Guidelines for cleaning and sterilization
  - Suitable also for high risk infection patients!
  - 3-pin rigid cranial fixation
  - For supine, prone, lateral or sitting position\*
  - Excellent stability and flexibility
- \* For sitting position Crossbar Adaptor is required, see additional accessories below

\* For sitting position Crossbar Adaptor is required

Includes the items below:

103 003 009	Skull Clamp Teflon
103 002 009	Doro Swivel Adaptor Teflon
103 001 009	Doro Adjustable Base Unit Teflon

**Pins available to the product**

