

EVOLUTION SUPINE POSITIONERS

Evolution Supine Positioners provides enhanced pressure relief with a proven anatomical design.

- Unique and innovative design that redistributes weight of the patient's body and provides enhanced pressure relief compared to regular Patient Positioners.
- The positioners are made from two components. The inside of the product is comfortable Visco-elastic foam. The outer coating is a seamless, non-porous Dual Coated Hi-Tech Membrane with high tear strength making it a very robust product.
- The outer coating has Excellent Abrasion Resistance and is Puncture/Burst Resistant.
- It has excellent flexibility that improves conformity to the skin.
- Provides Zero Interface Pressure under the Bony Prominences.
- High strength bond between foam and coating creates an ultra strong laminated structure for a long life product.
- The outer coating provides an excellent fluid-proof barrier without any seams or other bacteria harbouring features on the outer surface.
- Completely Latex free.
- Compatible with CT Scanners as these do not cause artefacts.
- Head Positioners provide excellent Neck Support and Head Stability compared to regular Head Rings.
- Foot Positioners provide excellent Foot Support with Zero Pressure under the Heel compared to regular Heel Pads.
- Competitively priced to make it a cost effective option against regular positioners.
- Two year product warranty.



EVOLUTION Supine Head Positioner

Raised neck area supporting cervical vertebrae and trapezoid muscle
Zero pressure on the Occipital area helps to prevent decubitus ulcers, pressure alopecia and fluid pooling
Fingered lateral support provide excellent pressure distribution and stable head support with comfortable fit



EVOLUTION Supine Foot Positioner

- Ergonomical support surface contoured to the back of the foot and the achilles tendon
- Raised sides protect the foot from lateral force and leg rotation
- Zero pressure on the calcaneus bone helps prevent decubitus ulcers developing

Smooth contoured surfaces allows maximum patient comfort
Non-slip base grips to operating table mattress
Air ventilation through the base maintains normothermia at the occipital bone
Concave base aids self adjusting motion and allows air ventilation
100066112 Medium Style 2
100066113 Large Style 2

- Smooth contoured surface allows maximum surface redistribution
 - Unique non-slip base grips to operating table mattress
 - Air ventilation through the base allows normothermia at the calcaneus bone
 - Concave base aids self adjustment motion and allows air ventilation
- 100066114 Medium

Product Code

100066112, 100066113, 100066114
