

# SPINE FRAME I 2805

Stable and convenient platform for patient positioning for spine surgeries such as laminectomies, decompressions and disc surgery with the Wilson Frame.

Includes frame, pads and a crank handle.

- Can be adjusted intra-operatively to the desired lordosis to get optimal surgical site exposure
- C-arm compatible
- Safe Patient weight 227 kg (500 lbs)



Bredd	460 mm
Längd	760 mm
Höjjustering	180 - 250 mm
Säker arbetsbelastning (SWL)	227 kg
Vikt	12,1 kg
Product Code	100012805, 12805